

See the future
SIUI

CTS-5500Plus

Portable Digital
B/W Ultrasound System




CTS-5500Plus

◎ Specification


- Monitor: 10-inch LCD monitor
- Depth: Scanning depth up to 252mm
- B Gain: 0~100dB
- Probe frequency range from 2.5 to 10.0 MHz
- Two probe connectors as standard
- USB port: 2
- Memory: 320G

◎ Feature


- Powerful digital beamforming technology
- Tissue Harmonic Imaging
- Unique high-definition zooming function
- IP one-key optimization
- Complete applications
- Comprehensive measurement package



Carotid



Hepatorenal section



MSK

SIUI

Shantou Institute of Ultrasonic Instruments Co., Ltd.

Add: #77, Jinsha Road, Shantou 515041, Guangdong, China

Tel: +86-754-88250150 Fax: +86-754-88251499

E-mail: siui@siui.com Website: <http://www.siui.com>

Specifications and appearance are subject to change without prior notice.

DCY2.782.EN.CTS-5500Plus_CY7A01