

## Programmable open system

- Real-time kinetic plots
- Levey-Jennings graphs for QC tracking
- Stores calibration and standard curves
- Stores and reports patient data

## Biochemistry Analyzer



3300 shown with optional keyboard and heat block

Stat Fax®, serving laboratory needs worldwide for more than 25 years

## Famous For Durability


- Accepts both round and square test tubes or flow cell
- Bichromatic optical system with 6 filters and space for 2 additional filters (340 nm - 700 nm)
- Calibration traceable to NIST standards
- Multiple calculation modes: factor, single standard, point-to-point, regression, and index calculation
- Reads, calculates, and prints results, 3 seconds per tube
- Store up to 120 tests in memory
- Stores over 500 patient results and 500 control results
- Optional incubation block at 37° and external keyboard

The compact design, large LCD screen and user-friendly interface make the Stat Fax® 3300 biochemistry analyzer an excellent choice for a variety of laboratory applications.

 **AWARENESS®**  
TECHNOLOGY, INC.  
[sales@awaretech.com](mailto:sales@awaretech.com)


# Stat Fax<sup>®</sup> 3300 Biochemistry Analyzer

## View kinetic plots in real time



See the curve, then calculate kinetic results based on the entire read time or the segment you choose.


## Stat Fax 3300 Stores calibration & standard curves



For each run, choose a new curve or use a stored curve with options to modify multi-point curves and revert to previous curve.

## Stat Fax 3300 Stores control values

With many QC options which assures linear kinetics and reliable results.



## Calculation Modes

Factor, Single Standard, Point-to-Point, Regression, Index Calculation



**Awareness Technology, Inc. USA & INTERNATIONAL:**  
 1935 S.W. Martin Hwy., Palm City, FL 34990, USA  
 tel.: +1 772 283 6540, fax: +1 772 283 8020  
 info@awaretech.com

**Awareness Technology, Inc. EUROPE:**  
 Franz-Siegel-Gasse 1, 2380, Perchtoldsdorf, Austria  
 tel.: +43 (1) 804 81 84, fax: +43 (1) 804 81 85  
 info@awaretech.eu

**Awareness Technology, Inc. ASIA & AFRICA**  
 Q3-146, SAIF Zone, 120543, Sharjah, United Arab Emirates  
 tel.: +971 (06) 557 80 58, fax: +971 (06) 557 80 59  
 info@awaretech.com



ISO 13485:2016

